

APD350Hu01 100µg
Active Interleukin 17 Receptor C (IL17RC)
Organism Species: Homo sapiens (Human)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

1st Edition (Apr, 2016)

[PROPERTIES]

Source: Prokaryotic expression.

Host: *E. coli*

Residues: Asn257~Asn477

Tags: N-terminal His-tag

Purity: >98%

Buffer Formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 0.05% sarcosyl and 5% trehalose.

Applications: Cell culture; Activity Assays; In vivo assays.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 5.8

Predicted Molecular Mass: 28.6kDa

Accurate Molecular Mass: 26kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

```
NHTQ QLPDCRGLEV WNSIPSCWAL PWLNVSADGD NVHLVLNVSE
EQHFGLSLYW NQVQGPPKPR WHKNLTGPQI ITLNHTDLVP CLCIQVWPLE
PDSVRTNICP FREDPRAHQN LWQAARLQLL TLQSWLLDAP CSLPAEAAALC
WRAPGGDPCQ PLVPPLSWEN VTVDKVLEFP LLKGHPNLCV QVNSSEKLQL
QECLWADSLG PLKDDVLLLE TRGPQDN
```

[ACTIVITY]

IL17RC (Interleukin-17 receptor C) is a receptor for IL17A and IL17F homodimers as part of a heterodimeric complex with IL17RA. This protein is expressed in nonhemopoietic tissues, and binds both IL-17A and IL-17F with similar affinities. A binding ELISA assay was conducted to detect the association of IL17RC with IL17RA. Briefly, IL17RC were diluted serially in PBS with 0.01%BSA (pH 7.4). Duplicate samples of 100uL were then transferred to IL17RA-coated microtiter wells and incubated for 2h at 37°C. Wells were washed with PBST and incubated for 1h with anti-IL17RC pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37°C. Finally, add 50µL stop solution to the wells and read at 450nm immediately. The binding activity of IL17RC with IL17RA was shown in Figure 1 and this effect was in a dose dependent manner.




Figure 1. The binding activity of IL17RC with IL17RA.

[IDENTIFICATION]




Figure 2. SDS-PAGE

Sample: Active recombinant IL17RC, Human




Figure 3. Western Blot

Sample: Recombinant IL17RC, Human;

Antibody: Rabbit Anti-Human IL17RC Ab (PAD350Hu01)