#### PAC130Ga01 Polyclonal Antibody to Noggin (NOG) **Organism Species: Chicken (Gallus)** Instruction manual

#### FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

[PRODUCT INFORMATION] Western Bolt kDa 75 Immunogen: NOG, Gallus 50 **Clonality:** Polyclonal 37 Host: Rabbit 25 ←24kDa Immunoglobulin Type: IgG 20 **Purification:** Affinity Chromatography. 15 Applications: WB, ICC, IHC-P, IHC-F, ELISA 10 Concentration: 200µg/mL **UOM:** 100µg Sample: Recombinant NOG, Gallus

### [IMMUNOGEN INFORMATION]

Immunogen: Recombinant NOG (Gln27~Cys223) expressed in E.coli. USCN Accession No.: RPC130Ga01 Sequence: The target protein is fused with N-terminal His-Tag and its sequence is listed below. MGHHHHHHSGSEF-QHYL HIRPAPSDNL PLVDLIEHPD PIFDPKEKDL NETLLRSLMG GHFDPNFMAM SLPEDRLGVD DLAELDLLLR QRPSGAMPGE IKGLEFYDGL QPGKKHRLSK KLRRKLQMWL WSQTFCPVLY TWNDLGSRFW PRYVKVGSCY SKRSCSVPEG MVCKPAKSVH LTILRWRCQR RGGQRCTWIP IQYPIIAECK CSC

9th Edition (Revised in Jul, 2013)

## [ANTIBODY SPECIFITY]

The antibody is a rabbit polyclonal antibody raised against NOG. It has been selected for its ability to recognize NOG in immunohistochemical staining and western blotting.

#### [APPLICATIONS]

Western blotting: 1:100-400 Immunocytochemistry in formalin fixed cells: 1:100-500 Immunohistochemistry in formalin fixed frozen section: 1:100-500 Immunohistochemistry in paraffin section: 1:50-200 Enzyme-linked Immunosorbent Assay: 1:100-200 Optimal working dilutions must be determined by end user.

# [<u>CONTENTS</u>]

**Form & Buffer:** Supplied as solution form in PBS, pH7.4, containing 0.02% NaN<sub>3</sub>, 50% glycerol.

### [ QUALITY CONTROL ]

**Content:** The quality control contains recombinant NOG (Gln27~Cys223) disposed in loading buffer.

Usage: 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.

5uL per well when used in enhanced chemilumescent (ECL). **Note:** The quality control is specifically manufactured as the positive control. Not used for other purposes.

Loading Buffer: 100mM Tris(pH8.8), 2% SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, NaN $_3$  0.02%.

# [<u>STORAGE</u>]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.

Designed by Cloud-Clone Corp., Assembled by Uscn Life Science Inc. ISO9001:2008; ISO13485:2003 11271 Richmond Avenue Suite H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Http://www.doud-clone.us.j.E-mail: mail@cloud-clone.us Export Processing Zone Building F, Wuhan, Hubel 430056, PRC | Toll free: 0086-800-880-0667 | Fax: 0086-27-8425-9551 | Http://www.uscnk.com | E-mail: mail@cloud-clone.us