

EPA363Ra61 100µg
Eukaryotic Platelet/Endothelial Cell Adhesion Molecule (PECAM1)
Organism Species: *Rattus norvegicus* (Rat)
Instruction manual

FOR RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[PROPERTIES]

Source: Eukaryotic expression

Host: 293F cell

Residues: Gln18~Lys589

Tags: N-terminal His Tag

Subcellular Location: Membrane

Purity: > 90%

Traits: Freeze-dried powder

Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 50µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 6.9

Predicted Molecular Mass: 66.2kDa

Accurate Molecular Mass: 120kDa as determined by SDS-PAGE reducing conditions.

Phenomenon explanation:

The possible reasons that the actual band size differs from the predicted are as follows:

1. Splice variants: Alternative splicing may create different sized proteins from the same gene.
2. Relative charge: The composition of amino acids may affect the charge of the protein.
3. Post-translational modification: Phosphorylation, glycosylation, methylation etc.
4. Post-translation cleavage: Many proteins are synthesized as pro-proteins, and then cleaved to give the active form.
5. Polymerization of the target protein: Dimerization, multimerization etc.

[USAGE]

Reconstitute in ddH₂O to a concentration of 0-0.5 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

```

                QEN SFTINSIHME SRPSWEVSNG QKLTQLCLVD
ISTTSKSRPQ HQVLFYKDDA LVYNVSSSEH TEFVIPQSR VFHAGKYKCT
VILNSKEKTT IEYQLTVNGV PMPEVTVDKK EVTEGGIVTV NCSMQEEKPP
IYFKIEKVEL GTKNVKLSRE KTSNMNFVLI EFPIDEEQDHL LVFRCQAGVL
SGIKMQTSEF IRSEYVTVQE FFSTPKFQIQ PPEMIEGNQ LHIKCSVQVA
HLAQEFPEII IQKDKAIVAT SKQSKEAVYS VMALVEHSGH YTCKVESNRI
SKASSILVNI TELFPRPKLE LSSRLDQGE MLDLSCSVSG APVANFTIQK
EETVLSQYQN FSKIAEERDS GLYSCTAGIG KVVKRSNLVP VQVCEMLSKP
RIFHDAKFEI IKGQIIGISC QSVNGTAPIT YRLLRAKSNF QTVQKNSNDP
VTFTDKPTRD MEYQCIVDNC HSHPEVRSEI LRVKVIAPVD EVTISILSGN
DVQSGDEMVL RCSVKEGTGP VTFQFYKEKE GRPFHEETVN DTQVFWHHEQ
TSKEQEGQYY CTAFNRASIV TSLRSGPLTV RVFLAPWKK
    
```

[IDENTIFICATION]




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.