

MAK477Mu21

Monoclonal Antibody to Transmembrane Protein 27 (TMEM27)

Organism Species: *Mus musculus* (Mouse)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Monoclonal antibody preparation

Host: Mouse

Antibody isotype: IgG1 Kappa

Purification: Protein A + Protein G affinity chromatography

Clone number: C7

Traits: Liquid

Concentration: 1mg/mL

UOM: 200μL

Cross Reactivity: Rat

Applications: WB; IHC

[IMMUNOGEN]

Immunogen: Recombinant TMEM27 (Ile12~Leu222) expressed in *E.coli*

Accession No.: RPK477Mu01

[APPLICATIONS]

Western blotting: 0.01-3μg/mL;

Immunohistochemistry: 5-20μg/mL;

Optimal working dilutions must be determined by end user.

[FORMULATION]

Form & Buffer: Supplied as solution form in 0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.



Store at 4°C for frequent use.

Aliquot and store at -20°C for 24 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no

obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[IDENTIFICATION]



DAB staining on IHC-P; Samples:	Western Blot; Sample: Rat Spleen
Mouse Kidney Tissue; Primary Ab:	lysate
20µg/ml Mouse Anti-Mouse TMEM27	Primary Ab: 3µg/ml Mouse Anti-Mouse
Antibody Second Ab: 2µg/mL HRP-	TMEM27 Antibody
Linked Caprine Anti-Mouse IgG	Second Ab: 0.2µg/mL HRP-Linked
Polyclonal Antibody (Catalog:	Caprine Anti-Mouse IgG Polyclonal
SAA544Mu19)	Antibody
	(Catalog: SAA544Mu19) Selected

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.