

**RPC517Hu01 200µg**

**Recombinant Hyaluronan Binding Protein 2 (HABP2)**

**Organism Species: *Homo sapiens (Human)***

***Instruction manual***

**FOR RESEARCH USE ONLY**

**NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES**

---

13th Edition (Revised in Aug, 2023)

## **[ PROPERTIES ]**

**Source:** Prokaryotic expression

**Host:** *E.coli*

**Residues:** Phe24~Phe560

**Tags:** N-terminal His Tag

**Subcellular Location:** Secreted

**Purity:** > 80%

**Traits:** Freeze-dried powder

**Buffer formulation:** PBS, pH7.4, containing 0.01% Sarcosyl, 5% Trehalose.

**Original Concentration:** 150µg/mL

**Applications:** Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

**Predicted isoelectric point:** 6.0

**Predicted Molecular Mass:** 63.9kDa

**Accurate Molecular Mass:** 75/55kDa as determined by SDS-PAGE reducing conditions.

### **Phenomenon explanation:**

The possible reasons that the actual band size differs from the predicted are as follows:

1. Splice variants: Alternative splicing may create different sized proteins from the same gene.
2. Relative charge: The composition of amino acids may affects the charge of the protein.
3. Post-translational modification: Phosphorylation, glycosylation, methylation etc.
4. Post-translation cleavage: Many proteins are synthesized as pro-proteins, and then cleaved to give the active form.
5. Polymerization of the target protein: Dimerization, multimerization etc.

## **[ USAGE ]**

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.


## **[ STORAGE AND STABILITY ]**

**Storage:** Avoid repeated freeze/thaw cycles.

Aliquot and store at -80°C for 12 months.

**[ SEQUENCE ]**

|             |            | FSLMSLL    | ESLDPDWTPD | QYDYSYEDYN |
|-------------|------------|------------|------------|------------|
| QEENTSSSTLT | HAENPDWYYT | EDQADPCQPN | PCEHGGDCLV | HGSTFTCSCL |
| APFSGNKCQK  | VQNTCKDNPC | GRGQCLITQS | PPYYRCVCKH | PYTGPSCSQV |
| VPVCRPNPCQ  | NGATCSRHKR | RSKFTCACPD | QFKGKFCEIG | SDDCYVGDGY |
| SYRGKMNRV   | NQHACLYWNS | HLLQENYNM  | FMEDAETHGI | GEHNFCRNP  |
| ADEKPWCFIK  | VTNDKVKWEY | CDVSACSAQD | VAYPEESPT  | PSTKLPGFDS |
| CGKTEIAERK  | IKRIYGGFKS | TAGKHPWQAS | LQSSLPLTIS | MPQGHFCGGA |
| LIHPCWVLTA  | AHCTDIKTRH | LKVVLGDQDL | KKEEFHEQSF | RVEKIFKYSH |
| YNERDEIPHN  | DIALCLKLPV | DGHCALESKY | VKTVCLPDGS | FPSGSECHIS |
| GWGVTETGKG  | SRQLLDAKVK | LIANTLCNSR | QLYDHMIDDS | MICAGNLQKP |
| GQDTCQGDSC  | GPLTCEKDGT | YYVYGIVSWG | LECGKRPVGY | TQVTKFLNWI |
| KATIKSESGF  |            |            |            |            |

[illegible]

**Figure . Gene Sequencing (extract)**




Figure. SDS-PAGE

**[ IMPORTANT NOTE ]**

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.