

RPE189Mu01 200µg

Recombinant Proprotein Convertase Subtilisin/Kexin Type 9 (PCSK9)

Organism Species: *Mus musculus* (Mouse)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Ser156~Gln694

Tags: N-terminal His and GST Tag

Subcellular Location: Cytoplasm

Purity: > 90%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% Sarcosyl, 5% Trehalose.

Original Concentration: 400µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 7.5

Predicted Molecular Mass: 87.3kDa

Accurate Molecular Mass: 87kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.


Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

TEEDRSPDGS	SQVEVYLLDT	SIQGAHREIE	GRVTITDFNS	VPEEDGTRFH	SIPWN LERIIPAWHQ
THLAGVVSGR	DAGVAKGTS	HSLRVLNCQG	KGTVSGTLIG	LEFIRKSQLI	RQASKCDSHG
PLAGGYSRIL	NAACRHLART	GVVLVAAAGN	FRDDACLYSP	ASAPEVITVG	QPSGPLVLL
LGTLTGNFGR	CVDLFAPGKD	IIGASSDCST	CFMSQSGTSQ	AAAHVAGIVA	ATNAQDQPV
LAELRQRLIH	FSTKDVINMA	WFPEDQQVLT	PNLVATLPPS	THETGGQLLC	RMLSREPTLT
TRTATATARC	APEEELLSCS	SFSRSGRRRG	DWIEAIGGQQ	VCKALNAFGG	EGVYAVARCC
LVPRANCSIH	NTPAARAGLE	THVHCHQKDH	VTGCSFHWE	VEDLSVRRQP	ALRSRRQPGQ
CVGHQAASVY	ASCCHAPGLE	CKIKEHGISG	PSEQVTVACE	AGWTLTGCVN	LPGASLTGGA
YSVDNLCVAR	VHDTARADRT	SGEATVAAAI	CCRSRPSAKA	SWVQ	

[IDENTIFICATION]



[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.