

RPC802Hu01 50µg

Recombinant Tensin 1 (TNS1)

Organism Species: *Homo sapiens (Human)*

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[**PROPERTIES**]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Ser4~Gly307

Tags: N-terminal His Tag

Subcellular Location: Cytoplasm

Purity: > 90%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% Sarcosyl, 5% Trehalose.

Original Concentration: 150µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 8.6

Predicted Molecular Mass: 38.3kDa

Accurate Molecular Mass: 38kDa as determined by SDS-PAGE reducing conditions.

[**USAGE**]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[**STORAGE AND STABILITY**]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[**SEQUENCE**]

SRTMEDS	CELDLVYVTE	RIIAVSFPST	ANEENFRSNL	REVAQMLKSK
HGGNYLLFNL	SERRPDITKL	HAKVLEFGWP	DLHTPALEKI	CSICKAMDTW
LNADPHNVVV	LHNKGNRGRI	GVVIAAYMHY	SNISASADQA	LDRFAMKRFY
EDKIVPIGQP	SQRRYVHYFS	GLLSGSIKMN	NKPLFLHHVI	MHGIPNFESK
GGCRPFLRIY	QAMQPVYTSG	IYNIPGDSQT	SVCITIEPGL	LLKGDILLKC
YHKKFRSPAR	DVIFRVQFHT	CAIHDLGUVF	GKEDLDDAFK	DDRFPEYGKV
EFVFSYG				

[IDENTIFICATION]





Figure . Gene Sequencing (extract)



[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.