

RPC036Hu01 50µg

Recombinant Transferrin (TF)

Organism Species: *Homo sapiens (Human)*

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Val361~Lys683

Tags: N-terminal His Tag

Subcellular Location: Secreted

Purity: > 80%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% Sarcosyl, 5% Trehalose.

Original Concentration: 500µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 5.1

Predicted Molecular Mass: 42.6kDa

Accurate Molecular Mass: 41kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

	VKWCALSHHE	RLKCEWSVN	SVGKIECVSA	ETTEDCIAKI
MNGEADAMSL	DGGFVYIAGK	CGLVPVLAEN	YNKSDNCEDT	PEAGYFAIAV
VKKSASDLTW	DNLKGKKSCH	TAVGRTAGWN	IPMGLLYNKI	NHCRFDEFFS
EGCAPGSKKD	SSLCKLCMGS	GLNLCEPNNK	EGYYGYTGAF	RCLVEKGDVA
FVKHQTPVQN	TGGKNPDPWA	KNLNEKDYEL	LCLDGTRKPV	EEYANCHLAR
APNHAVVTRK	DKEACVHKIL	RQQQHLFGSN	VTDCSGNFCL	FRSETKDLLF
RDDTVCLAKL	HDRNTYEKYL	GEEYVKAAGN	LRK	

[IDENTIFICATION]




Figure . Gene Sequencing (extract)




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.