

APB540Hu61 100µg
Active Cluster Of Differentiation 147 (CD147)
Organism Species: *Homo sapiens (Human)*
Instruction manual

FOR RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

1st Edition (Apr, 2016)

[PROPERTIES]

Source: Eukaryotic expression.

Host: 293F cell

Residues: Glu138~Ala323

Tags: N-terminal His-tag

Purity: >95%

Endotoxin Level: <1.0EU per 1µg (determined by the LAL method).

Buffer Formulation: PBS, pH7.4, containing 5% trehalose.

Applications: Cell culture; Activity Assays.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 4.8

Predicted Molecular Mass: 21.6kDa

Accurate Molecular Mass: 27-33kDa as determined by SDS-PAGE reducing conditions.

Phenomenon explanation:

The possible reasons that the actual band size differs from the predicted are as follows:

1. Splice variants: Alternative splicing may create different sized proteins from the same gene.
2. Relative charge: The composition of amino acids may affects the charge of the protein.
3. Post-translational modification: Phosphorylation, glycosylation, methylation etc.
4. Post-translation cleavage: Many proteins are synthesized as pro-proteins, and then cleaved to give the active form.
5. Polymerization of the target protein: Dimerization, multimerization etc.

[USAGE]

Reconstitute in 10mM PBS (pH7.6) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

```
EPG TVFTTVEDLG  
SKILLTCSLN DSATEVTGHR WLKGGVVLKE DALPGQKTEF KVDSDDQWGE  
YSCVFLPEPM GTANIQLHGP PRVKAVKSSE HINEGETAML VCKSESVPPP  
TDWAWYKITD SEDKALMNGS ESRFFVSSSQ GRSELHIENL NMEADPGQYR  
CNGTSSKGSD QAIITLRVRS HLA
```

[ACTIVITY]

Extracellular matrix metalloproteinase (MMP) inducer (EMMPRIN), also known as basigin and CD147, is a 44-66 kDa, variably N- and O-glycosylated, type I transmembrane protein that belongs to the immunoglobulin superfamily. Human EMMPRIN is 269 amino acids (aa) in length and contains a 24 aa signal sequence, a 183 aa extracellular domain (ECD), a 21 aa transmembrane (TM) segment and a 41 aa cytoplasmic tail. The ECD contains one C2-type and one V-type Ig-like domain. EMMPRIN is expressed in areas of tissue remodeling, including endometrium, placenta, skin, and regions undergoing angiogenesis. It is also expressed on cells with high metabolic activity, such as lymphoblasts, macrophages and particularly tumor cells.

A functional ELISA assay was conducted to detect the interaction of recombinant Human EMMPRIN/CD147 and recombinant rat Cyclophilin B (CYPB) . Briefly, biotin-linked CD147 were diluted serially in PBS, with 0.01% BSA (pH 7.4).

Duplicate samples of 100 μ l were then transferred to CYPB-coated microtiter wells and incubated for 1h at 37 $^{\circ}$ C . Wells were washed with PBST 3 times and incubation with Streptavidin-HRP for 30min, then wells were aspirated and washed 5 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37 $^{\circ}$ C . Finally, add 50 μ l stop solution to the wells and read at 450 nm immediately. The binding activity of CD147 and CYPB was shown in Figure 1, the EC50 for this effect is 53.34 μ g/mL.




Figure 1. The binding activity of CD147 and CYPB

[IDENTIFICATION]

CGAATTGGGATTCGGATATGGTACATGGTGAATGGTGAATGCTTAAAGCTGTTTGGTCAAGTGCAGGCTTAGTGCAGTGCAGGAGGAGTTCATGCTGGTGCATGCTGGGCTAGCTCCCTATAGTGAATGTTAAAGGCTTGGATCTTGAAGTCCGAGGCTGGATGGTCCGGTGTCTTATATGGAGGCTCAAMCAAGCGTGGATGG




Figure 2. Gene Sequencing (extract)




Figure 3. SDS-PAGE

Sample: Active recombinant CD147, Human

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.