APB743Hu01 50µg Active Signal Transducer And Activator Of Transcription 3 (STAT3) Organism Species: *Homo sapiens* (Human) *Instruction manual*

FOR RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

1st Edition (Apr, 2016)

[PROPERTIES]

Source: Prokaryotic expression. Host: *E. coli* Residues: Ile576~Pro678 Tags: N-terminal His-tag Purity: >95% Endotoxin Level: <1.0EU per 1µg (determined by the LAL method). Buffer Formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% sarcosyl, 5%Trehalose. Applications: Cell culture; Activity Assays. (May be suitable for use in other assays to be determined by the end user.) Predicted isoelectric point: 6.3 Predicted Molecular Mass: 12.9kDa Accurate Molecular Mass: 14kDa as determined by SDS-PAGE reducing conditions. [USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Cloud-Clone Corp.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

ILALW NEGYIMGFIS KERERAILST KPPGTFLLRF SESSKEGGVT FTWVEKDISG KTQIQSVEPY TKQQLNNMSF AEIIMGYKIM DATNILVSPL VYLYPDIP

[ACTIVITY]

Signal transducer and activator of transcription 3 (STAT3) is a member of the STAT protein family. STAT3 mediates the expression of a variety of genes in response to cell stimuli, and thus plays a key role in many cellular processes such as cell growth and apoptosis. STAT3 is phosphorylated by receptor-associated Janus kinases (JAK), form homo- or heterodimers, and translocate to the cell nucleus where they act as transcription activators. Besides, Epidermal Growth Factor Receptor 2 (EGFR2) has been identified as an interactor of STAT3, thus a binding ELISA assay was conducted to detect the interaction of recombinant human STAT3 and recombinant human EGFR2. Briefly, STAT3 were diluted serially in PBS, with 0.01% BSA (pH 7.4). Duplicate samples of 100µL were then transferred to EGFR2-coated microtiter wells and incubated for 2h at 37 °C . Wells were washed with PBST and incubated for 1h with anti-STAT3 pAb, then aspirated and washed 3 times. After incubation with HRP labelled secondary antibody, wells were aspirated and washed 3 times. With the addition of substrate solution, wells were incubated 15-25 minutes at 37 °C. Finally, add 50µL stop solution to the wells and read at 450nm immediately. The binding activity of STAT3 and EGFR2 was shown in Figure 1, and this effect was in a dose dependent manner.

Cloud-Clone Corp.




Figure 1. The binding activity of STAT3 with EGFR2

[IDENTIFICATION]




Figure 2. SDS-PAGE

Sample: Active recombinant STAT3, Human

Cloud-Clone Corp.




Figure 3. Western Blot Sample: Recombinant STAT3, Human; Antibody: Rabbit Anti-Human STAT3 Ab (PAB743Hu01)

[<u>IMPORTANT NOTE</u>]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.