


A91685Hu01
Polyclonal Antibody to Integrin Alpha M (ITGaM)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
 NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

[PRODUCT INFORMATION]

Immunogen: ITGaM
Clonality: Polyclonal
Host: Rabbit
Species Reactivity: Human
Ig type: Rabbit IgG
Purification: Antibodies are purified by target protein affinity chromatography.
Applications: WB, ICC, IHC-P, IHC-F
Form: Liquid
Size: 100µg



Sample: Recombinant human ITGaM

[IMMUNOGEN INFORMATION]

Immunogen: Recombinant human ITGaM (Asn301~Thr451) expressed in *E.coli*.
Molecular Weight: 43.7kDa
USCN accession No.: P91685Hu01

Sequence: The target protein is fused with two N-terminal Tags, His-tag and GST-tag and its sequence is listed below.

MRNKKFELGLEFPNLPYYIDGDVKLTQSMARIYIADKHNMLGGCPKERAIEISMLEGAVLDIRY
 GVSRIAYSKDFETLKVDLFLSKLPEMLKMFEDRLCHKTYLNGDHVTHPDFMLYDALDVVLYMDP
 MCLDAFPKLVCFKKRIEAIPIQIDKYLKSSKYIAWPLQGQWQATFGGGDHPPKSDGSTSGSGHHH
 HHHSAGLVPRGSTAIGMKETAAKFERQHMDSPDLGTLEVLFGGPLGSEF-NTIASKPPRD
 HVFQVNNFEA LKTIQNQLRE KIFAIEGTQT GSSSSFEHEM SQEGFSAAIT SNGPLLSTVG

Unique product Superb quality Client favorite Nicest service  **ISO9001:2008**;  **ISO13485:2003**; 

11271 Richmond Avenue, H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Email: mail@uscnk.us | Http://www.uscnk.us
 Export Processing Zone, Wuhan 430056, P.R.China | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-9551 | Email: mail@uscnk.com | Http://www.uscnk.com

SYDWAGGVFL YTSKEKSTFI NMTRVDSMDN DAYLGAAAAI ILRNRVQSLV LGAPRYQHIG
LVAMFRQNTG MWESNANVKG T

[ANTIBODY SPECIFICITY]

Anti ITGaM is a rabbit polyclonal antibody raised against human ITGaM. This antibody has been selected for its ability to recognize human ITGaM in immunohistochemical staining and western blotting, non cross-reactive with other members of the family.

[APPLICATIONS]

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Optimal working dilutions must be determined by end user.

[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH 7.4, containing 0.02%Na₃N, 50% glycerol.

[QUALITY CONTROL]

Content: The quality control contains recombinant human ITGaM (Asn301~Thr451) disposed in loading buffer.

Usage: 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.

5uL per well when used in enhanced chemiluminescent (ECL).

Note: The quality control is specifically manufactured as the positive control. Not used for other purposes.

Loading buffer: 100mM Tris(PH 8.8), 2%SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, Na₃N 0.02%.

[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.

Unique product Superb quality Client favorite Nicest service  ISO9001:2008;  ISO13485:2003; 

11271 Richmond Avenue, H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Email: mail@uscnk.us | Http://www.uscnk.us
Export Processing Zone, Wuhan 430056, P.R.China | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-9551 | Email: mail@uscnk.com | Http://www.uscnk.com