


A92538Mu01
Polyclonal Antibody to Islet Cell Autoantigen 1 (ICA1)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
 NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

[PRODUCT INFORMATION]

Immunogen: ICA1
Clonality: Polyclonal
Host: Rabbit
Species Reactivity: Mouse
Ig type: Rabbit IgG
Purification: Antibodies are purified by target protein affinity chromatography.
Applications: WB, ICC, IHC-P, IHC-F
Form: Liquid
Size: 100µg



Sample: Recombinant mouse ICA1

[IMMUNOGEN INFORMATION]

Immunogen: Recombinant mouse ICA1 (Met1~Gly257) expressed in *E. coli*.
Molecular Weight: 35.9kDa
USCN accession No.: P92538Mu01
Sequence: The target protein is fused with two N-terminal Tags, His-tag and S-tag and its sequence is listed below.
 MHHHHHHSSG LVPRGSGMKE TAAKFERQH MDSPDLGTDD DDKAMADIGS EF-
 MSGHKCYSWE LQDRFAQDKS VVNKMQQKYW ETKQAFIKAT GKKEDEHVVA
 SDADLDAKLE LFHSIQRCL DLSKAIVLYQ KRICFLSQEE NELGKFLRSQ GFQDKTRAGK
 MMQATGKALC FSSQQLALR NPLCRFHQEV ETRHRAISD TWLTVNRMEQ YRTEYRGALL

Unique product Superb quality Client favorite Nicest service  ISO9001:2008;  ISO13485:2003; 

WMKDV SQELD PDLYKQMEKF R KVQTQVRLA KKNFDK LKMD VCQKV DLLGA
SRCNLLSHML ATYQT TLLHF WEKTSHTMAA IHESFKG

[ANTIBODY SPECIFICITY]

Anti ICA1 is a rabbit polyclonal antibody raised against mouse ICA1. This antibody has been selected for its ability to recognize mouse ICA1 in immunohistochemical staining and western blotting, non cross-reactive with other members of the family.

[APPLICATIONS]

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Optimal working dilutions must be determined by end user.

[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH 7.4, containing 0.02% Na₃N, 50% glycerol.

[QUALITY CONTROL]

Content: The quality control contains recombinant mouse ICA1 (Met1~Gly257) disposed in loading buffer.

Usage: 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.

5uL per well when used in enhanced chemiluminescent (ECL).

Note: The quality control is specifically manufactured as the positive control. Not used for other purposes.

Loading buffer: 100mM Tris(PH 8.8), 2%SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, Na₃N 0.02%.

[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.

Unique product Superb quality Client favorite Nicest service  ISO9001:2008;  ISO13485:2003; 

11271 Richmond Avenue, H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Email: mail@uscnc.us | Http://www.uscnc.us
Export Processing Zone, Wuhan 430056, P.R.China | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-9551 | Email: mail@uscnc.com | Http://www.uscnc.com