

PAB480Ra01

Polyclonal Antibody to Mannose Binding Lectin (MBL)

Organism Species: Rattus norvegicus (Rat)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

[PRODUCT INFORMATION]

Immunogen: MBL, Rat

Clonality: Polyclonal

Host: Rabbit


Immunoglobulin Type: IgG

Purification: Affinity Chromatography.

Applications: WB, ICC, IHC-P, IHC-F, ELISA

Concentration: 200µg/mL

UOM: 100µg



Sample: Recombinant MBL, Rat

[IMMUNOGEN INFORMATION]

Immunogen: Recombinant MBL (Glu19~Asp244) expressed in *E.coli*.

Accession No.: RPB480Ra01

Sequence: The target protein is fused with N-terminal His-Tag and its sequence is listed below.

MGHHHHHHSGS-ET LTEGAQSSCP VIACSSPGLN GFPGKDGH DG AKGEKGEPGQ
GLRGLQGPPG KVGPA GPPGN PGSKGATGPK GDRGESVEFD TTNIDLEIAA LRSEL RAMRK
WVLLSMS ENV GK KYFMSSVR RMPLNRAKAL CSELQGT VAT PRNAEENRAI QNVAKDVAFL
GITDQRTENV FEDLTGN RVR YTNWNEGEPN NVGSGENCVV LLTNGKWN DV
PCSDSFLVVC EFSD

[ANTIBODY SPECIFICITY]

The antibody is a rabbit polyclonal antibody raised against MBL. It has been selected for its ability to recognize MBL in immunohistochemical staining and western blotting.

[APPLICATIONS]

Western blotting: 1:50-400

Immunocytochemistry in formalin fixed cells: 1:50-500

Immunohistochemistry in formalin fixed frozen section: 1:50-500

Immunohistochemistry in paraffin section: 1:10-100

Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.

[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN_3 , 50% glycerol.

[QUALITY CONTROL]

Content: The quality control is from rat liver extract known to contain the target protein MBL, disposed in loading buffer.

Usage: 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.
5uL per well when used in enhanced chemiluminescent (ECL).

Note: The quality control is specifically manufactured as the positive control. Not used for other purposes.

Loading Buffer: 100mM Tris(pH8.8), 2% SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, NaN_3 0.02%.

[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.