


A90705Hu01
Polyclonal Antibody to Mucin 2 (MUC2)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
 NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

[PRODUCT INFORMATION]

Immunogen: MUC2
Clonality: Polyclonal
Host: Rabbit
Species Reactivity: Human
Ig type: Rabbit IgG
Purification: Antibodies are purified by target protein affinity chromatography.
Applications: WB, ICC, IHC-P, IHC-F
Form: Liquid
Size: 100μg



Sample: Recombinant human MUC2

[IMMUNOGEN INFORMATION]

Immunogen: Recombinant human MUC2 (Tyr4480~Ser4690) expressed in *E. coli*.

Molecular Weight: 24.9kDa

USCN accession No.: P90705Hu01

Sequence: The target protein is fused with N-terminal His-Tag and its sequence is listed below.

MGHHHHHHSGSEF- Y CTGWGDPHYV TFDGLYYSYQ GNCTYVLVEE ISPSVDNFGV
 YIDNYHCDPN DKVSCPRTL I VRHETQEVLI KTVHMMPMQV QVQVNRQAVA LPYKKYGLEV
 YQSGINYVVD IPELGVLVS Y NGLSFSVRLP YHRFGNNTKG QCGTCTNTTS DDCILPSGEI
 VSNCEAAADQ WLVNDPSKPH CPHSSSTTKR PAVTVPGGGK TTPHKDCTPS

Unique product **Superb quality** **Client favorite** **Nicest service**  **ISO9001:2008**;  **ISO13485:2003**; 

[**ANTIBODY SPECIFICITY**]

Anti MUC2 is a rabbit polyclonal antibody raised against human MUC2. This antibody has been selected for its ability to recognize human MUC2 in immunohistochemical staining and western blotting, non cross-reactive with other members of the family.

[**APPLICATIONS**]

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Optimal working dilutions must be determined by end user.

[**CONTENTS**]

Form & Buffer: Supplied as solution form in PBS, pH 7.4, containing 0.02%Na₃N, 50% glycerol.

[**QUALITY CONTROL**]

Content: The quality control contains recombinant human MUC2 (Tyr4480~Ser4690) disposed in loading buffer.

Usage: 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.

5uL per well when used in enhanced chemilumescient (ECL).

Note: The quality control is specifically manufactured as the positive control. Not used for other purposes.

Loading buffer: 100mM Tris(PH 8.8), 2%SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, Na₃N 0.02%.

[**STORAGE**]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.