


**A96918Hu01**  
**Polyclonal Antibody to Uromodulin (UMOD)**  
***Instruction manual***

FOR IN VITRO USE AND RESEARCH USE ONLY  
 NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

**[ PRODUCT INFORMATION ]**

**Immunogen:** UMOD  
**Clonality:** Polyclonal  
**Host:** Rabbit  
**Species Reactivity:** Human  
**Ig type:** Rabbit IgG  
**Purification:** Antibodies are purified by target protein affinity chromatography.  
**Applications:** WB, ICC, IHC-P, IHC-F  
**Form:** Liquid  
**Size:** 100µg



*Sample: Recombinant human UMOD*

**[ IMMUNOGEN INFORMATION ]**

**Immunogen:** Recombinant human UMOD (Glu334~Ser589) expressed in *E.coli*.  
**Molecular Weight:** 29.8kDa  
**USCN accession No.:** P96918Hu01  
**Sequence:** The target protein is fused with N-terminal His-Tag and its sequence is listed below.

MGHHHHHHSGS-ECGANDM    KVS LGKCQLK    SLGFDKVFMY    LSDSRCSGFN  
 DRDNRDWVSV    VTPARDGPCG    TVLTRNETHA    TYSNTLYLAD    EIIIRD LNIK    INFACSYPLD  
 MKVSLKTALQ    PMVSALNIRV    GGTGMFTVRM    ALFQTPSYTQ    PYQGSSVTLS    TEAFLYVGTM

*Unique product Superb quality Client favorite Nicest service*  ISO9001:2008;  ISO13485:2003; 

LDGGDLRFAL LLMTNCRYATP SSNATDPLKY FIIQDRCPHT RDSTIQVVEN GESSQGRFSV  
QMFRFAGNYD LVYLHCEVYL CDTMNEKCKP TCSGTRFRS

### **[ ANTIBODY SPECIFICITY ]**

Anti UMOD is a rabbit polyclonal antibody raised against human UMOD. This antibody has been selected for its ability to recognize human UMOD in immunohistochemical staining and western blotting, non cross-reactive with other members of the family.

### **[ APPLICATIONS ]**

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Optimal working dilutions must be determined by end user.

### **[ CONTENTS ]**

**Form & Buffer:** Supplied as solution form in PBS, pH 7.4, containing 0.02% Na<sub>3</sub>N, 50% glycerol.

### **[ QUALITY CONTROL ]**

**Content:** The quality control contains recombinant human UMOD (Glu334~Ser589) disposed in loading buffer.

**Usage:** 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.

5uL per well when used in enhanced chemiluminescent (ECL).

Note: The quality control is specifically manufactured as the positive control. Not used for other purposes.

**Loading buffer:** 100mM Tris(PH 8.8), 2%SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, Na<sub>3</sub>N 0.02%.

### **[ STORAGE ]**

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.