

**A96894Hu01**

**Polyclonal Antibody to WNT1 Inducible Signaling Pathway Protein 2 (WISP2)**  
***Instruction manual***

FOR IN VITRO USE AND RESEARCH USE ONLY  
 NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

**[ PRODUCT INFORMATION ]**

**Immunogen:** WISP2

**Clonality:** Polyclonal

**Host:** Rabbit

**Species Reactivity:** Human


**Ig type:** Rabbit IgG

**Purification:** Antibodies are purified by target protein affinity chromatography.

**Applications:** WB, ICC, IHC-P, IHC-F

**Form:** Liquid

**Size:** 100µg



*Sample: Recombinant human WISP2*

**[ IMMUNOGEN INFORMATION ]**

**Immunogen:** Recombinant human WISP2 (Gln24~Phe250) expressed in *E.coli*.

**Molecular Weight:** 24.3kDa

**USCN accession No.:** P96894Hu01

**Sequence:** The target protein is fused with N-terminal His-Tag and its sequence is listed below.

MGHHHHHHSGSEF-QLCPTPC	TCPWPPPRCP	LGVPLVLDGC	GCCRVCARRL
GEPCDQLHVC	DASQGLVCQP	GAGPGGRGAL	CLLAEDDSSC
ETFQPHCSIR	CRCEDGGFTC	VPLCSEDVRL	PSWDCPHPRR
			VEVLGKCCPE

**Unique product Superb quality Client favorite Nicest service**  **ISO9001:2008**;  **ISO13485:2003**;  **CE**

11271 Richmond Avenue, H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Email: mail@uscnk.us | Http://www.uscnk.us  
 Export Processing Zone, Wuhan 430056, P.R.China | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-9551 | Email: mail@uscnk.com | Http://www.uscnk.com

WVCGQGGLG TQPLPAQGPQ FSGLVSSLPP GVPCPEWSTA WGPCSTTCGL  
GMATRVSNQN RFCRLETQRR LCLSRPCPPS RGRSPQNSAF

### **[ ANTIBODY SPECIFICITY ]**

Anti WISP2 is a rabbit polyclonal antibody raised against human WISP2. This antibody has been selected for its ability to recognize human WISP2 in immunohistochemical staining and western blotting, non cross-reactive with other members of the family.

### **[ APPLICATIONS ]**

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Optimal working dilutions must be determined by end user.

### **[ CONTENTS ]**

**Form & Buffer:** Supplied as solution form in PBS, pH 7.4, containing 0.02% Na<sub>3</sub>N, 50% glycerol.

### **[ QUALITY CONTROL ]**

**Content:** The quality control contains recombinant human WISP2 (Gln24~Phe250) disposed in loading buffer.

**Usage:** 10uL per well when 3,3'-Diaminobenzidine(DAB) as the substrate.  
5uL per well when used in enhanced chemiluminescent (ECL).

Note: The quality control is specifically manufactured as the positive control. Not used for other purposes.

**Loading buffer:** 100mM Tris(PH 8.8), 2%SDS, 200mM NaCl, 50% glycerol, BPB 0.01%, Na<sub>3</sub>N 0.02%.

### **[ STORAGE ]**

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.