Coud-Clone Corp.

EPB686Hu61 100ug Eukaryotic B-Cell Activating Factor (BAFF) Organism Species: Homo sapiens (Human) Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

11th Edition (Revised in May, 2016)

[PROPERTIES]

Source: Eukaryotic expression. Host: 293F cell Residues: Ala134~Leu285 Tags: Two Tags, His-tag and Fc-tag Homology: Mouse 60%, Rat 26% Tissue Specificity: Spleen, marrow. Subcellular Location: Cell membrane. Single-pass type II membrane protein. Secreted. **Purity:** >95% Endotoxin Level: <1.0EU per 1µg (determined by the LAL method). Traits: Freeze-dried powder Buffer Formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 5%Trehalose and Proclin300. Original Concentration: 200ug/mL Predicted isoelectric point: 4.9 Predicted Molecular Mass: 43.9kDa Accurate Molecular Mass: 44kDa as determined by SDS-PAGE reducing conditions. **Applications:** SDS-PAGE; WB; ELISA; IP; CoIP; EMSA; Reporter Assays; Purification; Amine Reactive Labeling. (May be suitable for use in other assays to be determined by the end user.)

Cloud-Clone Corp.

[<u>USAGE</u>]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

AVQGPEE TVTQDCLQLI ADSETPTIQK GSYTFVPWLL SFKRGSALEE KENKILVKET GYFFIYGQVL YTDKTYAMGH LIQRKKVHVF GDELSLVTLF RCIQNMPETL PNNSCYSAGI AKLEEGDELQ LAIPRENAQI SLDGDVTFFG ALKLL

[IDENTIFICATION]




Figure 1. Gene Sequencing (extract)

1304 Langham Creek Dr, Suite 226, Huustun, TX 77084, USA | 001-888-960-7402 | www.cloud-cloue.us | mail@cloud-cloue.as Export Processing Zone, Wuhan, Hubei 430056, PRC | 0086-000-880-0687 | www.cloud-cloue.com | mail@cloud-cloue.com

Coud-Clone Corp.

44
44
33
26
22
18
14
10

Figure 2. SDS-PAGE

1304 Langham Creek De, Suite 226, Houston, TX 77084, USA | 001-888-960-7402 | www.cloud-clone.ux | mailiz.cloud-clone.ax Export Processing Zone, Wuhan, Hubei 430056, PRC | 0086-300-880-0687 | www.cloud-clone.com | mailiz.cloud-clone.com