

RPJ516Mu02 50µg

Recombinant Lecithin Cholesterol Acyltransferase (LCAT)

Organism Species: Mus musculus (Mouse)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

[PROPERTIES]

Residues: Ile290~Pro432

Tags: N-terminal His-Tag

Accession: P16301

Host: *E. coli*

Subcellular Location: Secreted.

Purity: >95%

Endotoxin Level: <1.0EU per 1µg

(determined by the LAL method).


Formulation: Supplied as lyophilized form in PBS, pH7.4, containing 1mM DTT, 5% trehalose, 0.01% sarcosyl and preservative.

Predicted isoelectric point: 5.7

Predicted Molecular Mass: 17.9kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)



15% SDS-PAGE

[USAGE]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The sequence of the target protein is listed below.

I STPNFN YTVQ DFERFF TDLH FEEGWH MFLQ SRDLLER LPA PGVEVY CLYG VGRPTP HTYI
YDHNFP YKDP VAALYED GDD TVATRST ELG GQWQGR QSQP VHLLPM NETD
HLNMVFS NKT LEHINA ILLG AYRTPK SPAA SP

[REFERENCES]

1. Warden C.H., *et al.* (1989) J. Biol. Chem. 264:21573-21581.
2. Church D.M., *et al.* (2009) PLoS Biol. 7:E1000112-E1000112.
3. Subbaiah P.V., Liu M. (1996) J. Lipid Res. 37:113-122.
4. Ng D.S., *et al.* (2002) J. Biol. Chem. 277:11715-11720.