

**MSI291Rb11****Medium for Primary Rabbit Peripheral Sciatic Nerve Fibroblasts (PSNC)*****Instruction manual***

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

1st Edition (Revised in Sept, 2025)

**[ Description ]**

Medium for Primary Rabbit Peripheral Sciatic Nerve Fibroblasts, is a liquid medium free of mycoplasma, bacteria and fungi and contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals and 10% concentrations of fetal bovine serum . It is a complete medium designed for optimal growth of Primary Rabbit Peripheral Sciatic Nerve Fibroblasts *in vitro*.

**[ Components ]**

500 ml of Medium for Primary Rabbit Peripheral Sciatic Nerve Fibroblasts consists of basal medium, 50 ml of FBS, 5 ml of Fibroblast Growth Supplement and 5 ml of penicillin/streptomycin solution.

*Note: FBS, Fibroblast Growth Supplement and P/S are not pre-mixed in Primary Rabbit Peripheral Sciatic Nerve Fibroblasts Medium; they must be added separately to make the complete Primary Rabbit Peripheral Sciatic Nerve Fibroblasts Medium.*

**[ Storage ]**

Store the basal medium at 2-8°C, Fibroblast Growth Supplement, FBS and P/S solutions at -20°C.

Protect from light.

**[ Shipping ]**

Basal medium : Ice pack refrigerated transport. Fibroblast Growth Supplement, FBS and P/S solution: dry ice.

**[ Usage ]**

Fibroblast Growth Supplement, FBS and P/S solutions were thawed at room temperature, sprayed with 75% ethanol on bottles and tubes, and transferred to a sterile operating table. In a sterile field, Fibroblast Growth Supplement, FBS and P/S solutions were added to the basal medium and mixed evenly to obtain the recombinant complete medium, which was directly used for culture of Primary Rabbit Peripheral Sciatic Nerve Fibroblasts *in vitro* .

*Note: When stored in the dark at 2-8°C, the reconstituted complete medium is stable for three months.*

**[ Important note ]**

- In order to maintain the best use effect of this product, do not place it in room temperature or high temperature environment for a long time.
- This product is for scientific research use only. It is not for diagnostic, therapeutic, clinical, family and other purposes.