



MSI068Rb11

## Medium for Primary Rabbit Tracheal Smooth Muscle Cells (TSMC)

### *Instruction manual*

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

1st Edition (Revised in Feb, 2024)

### [ Description ]

Medium for Primary Rabbit Tracheal Smooth Muscle Cells, is a liquid medium free of mycoplasma, bacteria and fungi and contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals and low concentrations of fetal bovine serum (5%). It is a complete medium designed for optimal growth of Primary Rabbit Tracheal Smooth Muscle Cells *in vitro*.

### [ Components ]

500ml of Medium for Primary Rabbit Tracheal Smooth Muscle Cells consists of basal medium, 25ml of FBS, 5 ml Smooth Muscle Cell growth supplement and 5 ml of penicillin/streptomycin solution.

*Note: FBS, Smooth Muscle Cell growth supplement and P/S solution are not pre-mixed in Primary Rabbit Tracheal Smooth Muscle Cells Medium; they must be added separately to make the complete Primary Rabbit Tracheal Smooth Muscle Cells Medium.*

### [ Storage ]

Store the basal medium at 4°C, FBS, Smooth Muscle Cell growth supplement and P/S solutions at -20°C.

Protect from light.

### [ Shipping ]

Basal medium: room temperature. FBS, Smooth Muscle Cell growth supplement and P/S solutions: dry ice.

### [ Usage ]

Smooth Muscle Cell growth supplement, FBS and P/S solutions were thawed at room temperature, sprayed with 75% ethanol on bottles and tubes, and transferred to a sterile operating table. In a sterile field, Smooth Muscle Cell growth supplement, FBS and P/S solutions were added to the basal medium and mixed evenly to obtain the recombinant complete medium, which was directly used for culture of Primary Rabbit Tracheal Smooth Muscle Cells *in vitro*.

*Note: When stored in the dark at 4°C, the reconstituted complete medium is stable for three months.*

### [ Important note ]

- In order to maintain the best use effect of this product, do not place it in room temperature or high temperature environment for a long time.
- This product is for scientific research use only. It is not for diagnostic, therapeutic, clinical, family and other purposes.