

MAA310Bo21

Monoclonal Antibody to Apolipoprotein H (APOH)

Organism Species: *Bos taurus*; Bovine (Cattle)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

[**PRODUCT INFORMATION**]

Immunogen: APOH, Bovine

Clonality: Monoclonal

Host: Mouse


Immunoglobulin Type: IgG

Purification: Affinity Chromatography.

Applications: WB, ICC, IHC-P, IHC-F, ELISA

Concentration: 500µg/mL

UOM: 200µg



Sample: Recombinant APOH, Bovine

[**IMMUNOGEN INFORMATION**]

Immunogen: Recombinant APOH (Arg21~Cys345) expressed in *E.coli*.

USCN Accession No.: RPA310Bo01

Sequence: The target protein is fused with N-terminal His-Tag and its sequence is listed below.

MGHHHHHHS~~G~~SEF-RTCPKPDEL~~P~~ FSTVVPLKRT YEPGEQIVFS CQPGYVSRGG
IRRFTCP~~L~~TG LWPINTLKCM PRVCPFAGIL EN~~G~~TVRYTTF EYPNTISFSC HTGFYLKGAS
SAKCTEEGKW SPDLPVCAPI TCPPPIPKF ASLSVYKPLA GNNSFYGSKA VFKCLPHHAM
FGNDTVTCTE HGNWTQLPEC REVRCPFPSR PDNGFVNHPA NPVLYYKDTA
TFGCHETYSL DGPEEVECSK FGNWSAQPSC KASCKLSIKR ATVIYEGERV AIQNKFKNGM
LHGQKVSFFC KHKEKKCSYT EDAQCIDGTI EIPKCFKEHS SLAFWKTDAS DVKPC

[ANTIBODY SPECIFICITY]

The antibody is a mouse monoclonal antibody raised against APOH. It has been selected for its ability to recognize APOH in immunohistochemical staining and western blotting.

[APPLICATIONS]

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.

[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.

[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.