

**MAB365Ra22**

**Monoclonal Antibody to Glucokinase Regulatory Protein (GKRP)**

**Organism Species: Rattus norvegicus (Rat)**

***Instruction manual***

FOR IN VITRO USE AND RESEARCH USE ONLY  
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

**[ PRODUCT INFORMATION ]**

**Immunogen:** GKRP, Rat

**Clonality:** Monoclonal

**Clone number:** F8

**Host:** Mouse


**Immunoglobulin Type:** IgG

**Purification:** Affinity Chromatography.

**Applications:** WB, ICC, IHC-P, IHC-F, ELISA

**Concentration:** 500µg/mL

**UOM:** 200µg



*Sample: Recombinant GKRP, Rat*

**[ IMMUNOGEN INFORMATION ]**

**Immunogen:** Recombinant GKRP (Val90~Val286) expressed in *E.coli*.

**Accession No.:** RPB365Ra01

**Sequence:** The target protein is fused with N-terminal His-Tag and its sequence is listed below.

MGHHHHHSGSEF-V QEVLKEPDGG LVVLSGGGTS GRMAFLMSVS FNQLMKGLGQ  
KPLYTYLIAG GDRSVVASRE QTEDSALHGI EELKKVAAGK KRVVIGISV GLSAPFVAGQ  
MDYCMDNTAV FLPVLVGFNP VSMARNDPIE DWRSTFRQVA ERMQKMQEKQ  
EAFVLNPAIG PEGLSGSSRM KGGGATKILL ETLLAAHKT VDQGVV

## **[ ANTIBODY SPECIFICITY ]**

The antibody is a mouse monoclonal antibody raised against GGRP. It has been selected for its ability to recognize GGRP in immunohistochemical staining and western blotting.

## **[ APPLICATIONS ]**

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.

## **[ CONTENTS ]**

**Form & Buffer:** Supplied as solution form in PBS, pH7.4, containing 0.02%  $\text{NaN}_3$ , 50% glycerol.

## **[ STORAGE ]**

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.