

MAC007Mu21

Monoclonal Antibody to Interleukin 34 (IL34)

Organism Species: *Mus musculus* (Mouse)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

[PRODUCT INFORMATION]

Immunogen: IL34, Mouse

Clonality: Monoclonal

Clone number: A10

Host: Rat


Immunoglobulin Type: IgG

Purification: Affinity Chromatography.

Applications: WB, ICC, IHC-P, IHC-F, ELISA

Concentration: 500µg/mL

UOM: 200µg



Sample: Recombinant IL34, Mouse

[IMMUNOGEN INFORMATION]

Immunogen: Recombinant IL34 (Asn21~Pro235) expressed in *E. coli*.

Accession No.: RPC007Mu01

Sequence: The target protein is fused with N-terminal His-Tag and its sequence is listed below.

MGHHHHHHSGSEF-NENLEIWTLT QDKECDLTGY LRGKLQYKNR LQYMKHYFPI
NYRIAVPYEG VLRVANITRL QKAHVSRERL RYLWVLVSLN ATEVMDVLL EGHPVSWKYLQ
EVQTLLENVQ RSLMDVEIGP HVEAVLSLLS TPGLSLKLVK PKALLDNCFR VMELLYCSCC
KQSPILKWQD CELPRLPHPS PGSLMQCTAT NVYPLSRQTP TSLPGSPSSS HGSLP

[ANTIBODY SPECIFICITY]

The antibody is a rat monoclonal antibody raised against IL34. It has been selected for its ability to recognize IL34 in immunohistochemical staining and western blotting.

[APPLICATIONS]

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.

[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.

[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.