

**MAB760Hu22**

**Monoclonal Antibody to Regenerating Islet Derived Protein 1 Alpha (REG1a)**

**Organism Species: Homo sapiens (Human)**

***Instruction manual***

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

## **[ PRODUCT INFORMATION ]**

**Immunogen:** REG1a, Human

**Clonality:** Monoclonal

**Clone number:** B3

**Host:** Mouse


**Immunoglobulin Type:** IgG

**Purification:** Affinity Chromatography.

**Applications:** WB, ICC, IHC-P, IHC-F, ELISA

**Concentration:** 500µg/mL

**UOM:** 200µg



*Sample: Recombinant REG1a, Human*

## **[ IMMUNOGEN INFORMATION ]**

**Immunogen:** Recombinant REG1a (Arg33~Asn166) expressed in *E.coli*.

**Accession No.:** RPB760Hu01

**Sequence:** The target protein is fused with N-terminal His-Tag and its sequence is listed below.

MGHHHHHHSGS-RISCP EGT NAYRSYCYYF NEDRETWVDA DLYCQNMNSG  
NLVSVLTQAE GAFVASLIKE SGTDDFNWVI GLHDPKKNRR WHWSSGSLVS YKSWGIGAPS  
SVNPGYCVSL TSSTGFQKWK DVPCEDKFSF VCKFKN

## **[ ANTIBODY SPECIFICITY ]**

The antibody is a mouse monoclonal antibody raised against REG1a. It has been selected for its ability to recognize REG1a in immunohistochemical staining and western blotting.

## **[ APPLICATIONS ]**

Western blotting: 1:100-400

Immunocytochemistry in formalin fixed cells: 1:100-500

Immunohistochemistry in formalin fixed frozen section: 1:100-500

Immunohistochemistry in paraffin section: 1:50-200

Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.

## **[ CONTENTS ]**

**Form & Buffer:** Supplied as solution form in PBS, pH7.4, containing 0.02% NaN<sub>3</sub>, 50% glycerol.

## **[ STORAGE ]**

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.