CPG918Hu21 100µg OVA Conjugated Uromodulin (UMOD) Organism Species: Homo sapiens (Human) *Instruction manual*

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

Coud-Clone Corp.

[PROPERTIES]

Source: Protein Conjugation Original Peptide Residues: Asp363-Val378 Purity: >90% Traits: Freeze-dried powder. Buffer formulation: PBS, pH7.4. Applications: Immunogen; Coating Antigen; Blocking Peptide; ELISA.

[<u>USAGE</u>]

Reconstitute in PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[<u>SEQUENCE</u>]

DSRCSGFN DRDNRDWV

Coud-Clone Corp.

[IDENTIFICATION]




Figure 1. HPLC of the Original Peptide




Figure 2. MASS Spectrometry of the Original Peptide

1304 Langham Creek Dr. Soite 226, Houston, TX 77084, USA | 001-888-960-7402 | www.cloud-chote.os | mail@cloud-clote.os Export Processing Zone, Wuhan, Hubei 430056, PRC | 0086-008-96087 | www.cloud-clote.com | mail@cloud-clote.com