

PAA804Mi01

Polyclonal Antibody to Glucagon Like Peptide 1 (GLP1)

Organism Species: Multi-species

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

[PRODUCT INFORMATION]

Immunogen: GLP1

Clonality: Polyclonal

Host: Rabbit

Immunoglobulin Type: IgG

Purification: Affinity Chromatography.

Applications: WB, ICC, IHC-P, IHC-F, ELISA

Concentration: 200µg/mL

UOM: 200µg

[IMMUNOGEN INFORMATION]

Immunogen: Synthetic Peptide GLP1.

Accession No.: CPA804Mi21

Sequence: The target peptide,s sequence is listed below.

HAEGTFTSDVSSYLEGQAAKEFIAWLVKGRG

[ANTIBODY SPECIFICITY]

The antibody is a rabbit polyclonal antibody raised against GLP1. It has been selected for its ability to recognize GLP1 in immunohistochemical staining and western blotting.

Note: Full length human GLP1 is identical in sequence to mouse, rat, porcine, bovine and canine. 100% cross-reactivity of GLP1 was observed among human, mouse, rat, porcine, bovine and canine.

[APPLICATIONS]

Western blotting: 1:50-400

Immunocytochemistry in formalin fixed cells: 1:50-500

Immunohistochemistry in formalin fixed frozen section: 1:50-500

Immunohistochemistry in paraffin section: 1:10-100

Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.

[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN₃, 50% glycerol.

[STORAGE]

Store at 4°C for frequent use. Stored at -20°C to -80°C in a manual defrost freezer for one year without detectable loss of activity. Avoid repeated freeze-thaw cycles.