

RPC992Mu01 200µg

Recombinant Glutathione Peroxidase 3, Plasma (GPX3)

Organism Species: *Mus musculus* (Mouse)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[**PROPERTIES**]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Arg3~Tyr72

Tags: N-terminal His Tag

Subcellular Location: Secreted

Purity: > 95%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 400µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 8.7

Predicted Molecular Mass: 11.4kDa

Accurate Molecular Mass: 12kDa as determined by SDS-PAGE reducing conditions.

[**USAGE**]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[**STORAGE AND STABILITY**]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[**SEQUENCE**]

RILRASCL LSLLLAGFVP PGRGQEKSKT DCHGGMSGTI YEYGALTIDG
EEYIPFKQYA GKYILFVNVA SY

[IDENTIFICATION]

XCGGATCC TCCGGGGCATCC TGCC TTCGTGTCCTGC TCC TGGCCGGGT TTTGTTCCGCGGGCCGGGG #CAAGAGAGTCTAAGACAGACTGCCATGGCGGTATGAGTGCTACCATCTACGAGTATGGAGCCCTC#CCATCGATGGGGAGG
 R I L R A S C L L S L L L A G F V P P G R G Q E K S K T D C H G G M S G T I Y E Y G A L T I D G E E





Figure. Gene Sequencing (Extract)



[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.