

**RPA822Mu01 1mg**

**Recombinant Heat Shock Protein 60 (Hsp60)**

**Organism Species: *Mus musculus* (Mouse)**

***Instruction manual***

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

---

13th Edition (Revised in Aug, 2023)

## **[ PROPERTIES ]**

**Source:** Prokaryotic expression

**Host:** *E.coli*

**Residues:** Lys28~Phe573

**Tags:** N-terminal His Tag

**Subcellular Location:** Mitochondrion

**Purity:** > 80%

**Traits:** Freeze-dried powder

**Buffer formulation:** 20mM Tris, 150mM NaCl, pH8.0, containing 0.01% skl, 5%Trehalose.

**Original Concentration:** 150µg/mL

**Applications:** Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

**Predicted isoelectric point:** 5.9

**Predicted Molecular Mass:** 61.4kDa

**Accurate Molecular Mass:** 65kDa as determined by SDS-PAGE reducing conditions.

### **Phenomenon explanation:**

The possible reasons that the actual band size differs from the predicted are as follows:

- 1.Splice variants: Alternative splicing may create different sized proteins from the same gene.
2. Relative charge: The composition of amino acids may affects the charge of the protein.
3. Post-translational modification: Phosphorylation, glycosylation, methylation etc.
4. Post-translation cleavage: Many proteins are synthesized as pro-proteins, and then cleaved to give the active form.
5. Polymerization of the target protein: Dimerization, multimerization etc.

## **[ USAGE ]**

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

## **[ STORAGE AND STABILITY ]**

**Storage:** Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

### [ SEQUENCE ]

```

                                KDV KFGADARALM LQGVDLLADA
VAVTMGPKGR TVIIEQSWGS PKVTKDGVTV AKSIDLKDKY KNIGAKLVQD
VANNTNEEAG DGTITATVLA RSIAGEGFEK ISKGANPVEI RRGVMLAVDA
VIAELKKQSK PVTTPPEIAQ VATISANGDK DIGNIISDAM KKVGRKGVIT
VKDGKTLNDE LEIIEGMKFD RGYISPYFIN TSKGQKCEFQ DAYVLLSEKK
ISSVQSIVPA LEIANHRKP LVIIAEDVDG EALSTLVLR LKVGQLQVAV
KAPGFGDNRK NQLKDAMIAT GGAVFGEEGL NLNLEDVQAH DLGKVGEVIV
TKDDAMLLKG KGDKAHIEKR IQEITEQLDI TTSEYEKEKL NERLAKLSDG
VAVLKVGGS DVEVNEKKDR VTDALNATRA AVEEGIVLGG GCALLRCIPA
LDSLKPANED QKIGIEIIKR ALKIPAMTIA KNAGVEGSLI VEKILQSSSE
VGYDAMLGDF VNMVEKGIID PTKVVRTALL DAAGVASLLT TAEAVVTEIP
KEEKDPGMGA MGGMGGGMGG GMF

```

### [ IDENTIFICATION ]




Figure. SDS-PAGE

### [ IMPORTANT NOTE ]



The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.