

RPB995Hu61 100µg

Recombinant Hemojuvelin (HJV)

Organism Species: Homo sapiens (Human)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

[PROPERTIES]

Residues: Gln36~Asp400

Tags: Two Tags, His-tag and Fc-tag

Accession: Q6ZVN8

Host: 293F cell

Subcellular Location: Cell membrane;

Lipid-anchor, GPI-anchor.

Purity: >95%

Endotoxin Level: <1.0EU per 1µg

(determined by the LAL method).


Formulation: Supplied as lyophilized form in PBS, pH7.4, containing 5% trehalose, 0.01% sarcosyl.

Predicted isoelectric point: 7.7

Predicted Molecular Mass: 66.1kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)



[USAGE]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[**STORAGE AND STABILITY**]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[**SEQUENCES**]

The sequence of the target protein is listed below.

QCKIL RCNAEYVSST LSLRGGGSSG ALRGGGGGGR GGGVGSGGLC RALRSYALCT
RRTARTCRGD LAFHSAVHGI EDLMIQHNCs RQGPTAPPPP RGPALPGAGS GLPAPDPCDY
EGRFSRLHGR PPGFLHCASF GDPHVRSFHH HFHTCRVQGA WPLLDNDFLF
VQATSSPMAL GANATATRK LTIIFKNMQEC IDQKVYQAEV DNL PVAFEDG SINGGDRPGG
SSLSIQTANP GNHVEIQAAY IGTTIIIRQT AGQLSFSIKV AEDVAMAFSA EQDLQLCVGG
CPPSQRLSRS ERNRRGAI TI DTARRLCKEG LPVEDAYFHS CVFDVLISGD PNFTVAAQAA
LEDARAFLPD LEKLHLFPSD