

**RPC517Hu01 50µg**

**Recombinant Hyaluronan Binding Protein 2 (HABP2)**

**Organism Species: *Homo sapiens (Human)***

***Instruction manual***

**FOR RESEARCH USE ONLY**

**NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES**

---

13th Edition (Revised in Aug, 2023)

## **[ PROPERTIES ]**

**Source:** Prokaryotic expression

**Host:** *E.coli*

**Residues:** Phe24~Phe560

**Tags:** N-terminal His Tag

**Subcellular Location:** Secreted

**Purity:** > 90%

**Traits:** Freeze-dried powder

**Buffer formulation:** 100mMNaHCO<sub>3</sub>, 500mMNaCl, pH8.3, containing 1mM EDTA, 1mM DTT, 0.01% Sarcosyl, 5% Trehalose and Proclin300.

**Original Concentration:** 200µg/mL

**Applications:** Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

**Predicted isoelectric point:** 6.0

**Predicted Molecular Mass:** 63.9kDa

**Accurate Molecular Mass:** 70kDa as determined by SDS-PAGE reducing conditions.

### **Phenomenon explanation:**

The possible reasons that the actual band size differs from the predicted are as follows:

1. Splice variants: Alternative splicing may create different sized proteins from the same gene.
2. Relative charge: The composition of amino acids may affects the charge of the protein.
3. Post-translational modification: Phosphorylation, glycosylation, methylation etc.
4. Post-translation cleavage: Many proteins are synthesized as pro-proteins, and then cleaved to give the active form.
5. Polymerization of the target protein: Dimerization, multimerization etc.

## **[ USAGE ]**

Reconstitute in ddH<sub>2</sub>O to a concentration of 0.1-1.0 mg/mL. Do not vortex.

## **[ STORAGE AND STABILITY ]**

**Storage:** Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.


### [ SEQUENCE ]

```

                                FSLMSLL ESLDPDWTPD QYDYSYEDYN
QEENTSSTLT HAENPDWYYT EDQADPCQPN PCEHGGDCLV HGSTFTCSCL
APFSGNKCQK VQNTCKDNPC GRGQCLITQS PPYYRCVCKH PYTGPSCSQV
VPVCRPNPCQ NGATCSRHKR RSKFTCACPD QFKGKFCEIG SDDCYVGDGY
SYRGKMNRTV NQHACLYWNS HLLLQENYNM FMEDAETHGI GEHNFCRNPD
ADEKPWCFIK VTNDKVKWEY CDVSACSAQD VAYPEESPT E PSTKLPGFDS
CGKTEIAERK IKRIYGGFKS TAGKHPWQAS LQSSLPLTIS MPQGHFCGGA
LIHPCWVLTA AHCTDIKTRH LKVVLGDQDL KKEEFHEQSF RVEKIFKYSH
YNERDEIPHN DIALCLKPV DGHCALESKY VKTVCLPDGS FPSGSECHIS
GWGVTETGKG SRQLLDKVK LIA NTLCNSR QLYDHMIDDS MICAGNLQKP
GQDTCQGDSG GPLTCEKDGT YYVYGIVSWG LECGKRPVY TQVTKFLNWI
KATIKSESGF

```

### [ IDENTIFICATION ]



**Figure . Gene Sequencing (extract)**




Figure. SDS-PAGE

### [ **IMPORTANT NOTE** ]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.