

RPA108Hu61 100µg

Recombinant Osteoprotegerin (OPG)

Organism Species: Homo sapiens (Human)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

[PROPERTIES]

Residues: Tyr31~Leu401

Tags: A his-tag

Accession: O00300

Host: 293F cell

Subcellular Location: Secreted.

Purity: >95%

Endotoxin Level: <1.0EU per 1µg

(determined by the LAL method).


Formulation: Supplied as lyophilized form in PBS, pH7.4, containing 5% trehalose, 0.01% sarcosyl.

Predicted isoelectric point: 8.7

Predicted Molecular Mass: 44.1kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)



15% SDS-PAGE

[USAGE]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The sequence of the target protein is listed below.

YDEETSHQLL CDKCPPGYL KQHCTAKWKT VCAPCPDHYY TDSWHTSDEC LYCSPVCKEL
QYVKQECNRT HNRVCECKEG RYLEIEFCLK HRSCPPGFGV VQAGTPERN
VCKRCPDGFF SNETSSKAPC RKHTNCSVFG LLLTQKGNAT HDNICS GNSE STQKCGIDVT
LCEEAFFRFA VPTKFTPNWL SVLVDNLPQT KVNAESVERI KRQHSSQEQT FQLLKLWKHQ
NKDQDIVKKI IQDIDLCENS VQRHIGHANL TFEQLRSLME SLPGKKVGAE DIEKTIKACK
PSDQILKLLS LWRIKNGDQD TLKGLMHALK HSKTYHFPKT VTQSLKKTIR FLHSFTMYKL
YQKL FLEMIG NQVQSVKISC L