

RPC919Hu02 50µg

Recombinant Platelet Derived Growth Factor D (PDGFD)

Organism Species: *Homo sapiens (Human)*

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Arg19~Arg370

Tags: N-terminal His Tag

Subcellular Location: Secreted

Purity: > 90%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 8.4

Predicted Molecular Mass: 44.4kDa

Accurate Molecular Mass: 45kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

```

RD TSATPQSASI KALRNANLRR DESNHLTDLY
RRDETIQVKG NGYVQSPRFP NSYPRNLLLW WRLHSQENTR IQLVFDNQFG
LEEAENDICR YDFVEVEDIS ETSTIIRGRW CGHKEVPPRI KSRTNVIKIT
FKSDDYFVAK PGFKIYYSLL EDFQAAAASE TNWESVTSSI SGVSYNPSV
TDPTLIADAL DKKIAEFDTV EDLLKYFNPE SWQEDLENNY LDTPRYRGRS
YHDRKSKVDL DRLNDDAKRY SCTPRNYSVN IREELKLANV VFFPRCLLVQ
RCGGNCGCGT VNWRSCTCNS GKTVKKYHEV LQFEPGHIKR RGRAKTMALV
DIQLDHERC DCICSSRPPR
  
```

[IDENTIFICATION]




Figure . Gene Sequencing (extract)




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.