

RPM090Hu01 50µg

Recombinant Ribosomal Protein S6 Kinase Alpha 5 (RPS6Ka5)

Organism Species: *Homo sapiens (Human)*

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Glu2~Phe348

Tags: N-terminal His Tag

Tissue Specificity: Heart, Brain, Placenta

Subcellular Location: Nucleus, Cytoplasm

Purity: > 95%

Traits: Freeze-dried powder

Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 5.9

Predicted Molecular Mass: 42.5kDa

Accurate Molecular Mass: 44kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the

expiration date under appropriate storage condition.

[SEQUENCE]

EEEGSSGG AAGTSADGGD GGEQLLTVKH ELRTANLTGH AEKVGIENFE
LLKVLGTGAY GKVFLVRKIS GHDTGKLYAM KVLKKATIVQ KAKTTEHTRT
ERQVLEHIRQ SPFLVTLHYA FQTETKLHLI LDYINGGELF THLSQRERFT
EHEVQIYVGE IVLALEHLHK LGIIYRDIKL ENILLDSNGH VVLTDFGLSK
EFVADETERA YSFCGTIEYM APDIVRGGDS GHDKAVDWWWS LGVLMYELLT
GASPFTVDGE KNSQAEISRR ILKSEPPYPQ EMSALAKDLI QRLLMKDPKK
RLGCGPRDAD EIKEHLFFQK INWDDLAACK VPAPFKPVIR DELDVSNF

[IDENTIFICATION]




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for in vitro and research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.