

P90367Mu01 100μg
Retinol Binding Protein 3, Interstitial (RBP3)
Organism: Mus musculus (Mouse)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
 NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

6th Edition (Revised in March, 2013)

[PROPERTIES]

Residues: Gly18~Leu320 (Accession # P49194),
 with N-terminal His-Tag.

Host: *E. coli*

Subcellular Location: Secreted, extracellular
 space, extracellular matrix, interphotoreceptor
 matrix.

Purity: >95%

Endotoxin Level: <1.0EU per 1μg
 (determined by the LAL method).


Formulation: Supplied as lyophilized form in PBS,
 pH7.4, containing 5% sucrose, 0.01% sarcosyl.

Predicted isoelectric point: 5.4

Predicted Molecular Mass: 34.5kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)



[USAGE]

Reconstitute in sterile PBS, pH7.2-pH7.4.

Unique product *Superb quality* *Client favorite* *Nicest service*  ISO9001:2008;  ISO13485:2003; 

11271 Richmond Avenue, H104, Houston, TX 77082, USA | Toll free: 001-888-960-7402 | Fax: 001-832-538-0088 | Email: mail@uscnk.us | Http://www.uscnk.us
 Export Processing Zone, Wuhan 430056, P.R.China | Toll free: 0086-800-880-0687 | Fax: 0086-27-8425-9551 | Email: mail@uscnk.com | Http://www.uscnk.com

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with N-terminal His-Tag, its sequence is listed below.

MGHHHHHHSGSEF-GPT HLFQPSLVLD MAKILLDNYC FPENLMGMQA AIEQAMKSHE
ILGISDPQTL AQVLTAGVQS SLSDPRLFIS YEPSTLEAPQ QAPVLTNLTR EELLAQIQRN
IRHEVLEGNV GYLRVDDLPG QEVLSSELGEF LVSHVWRQLM GTSSLVLDLR HCSSGGHFSGI
PYVISYLHPG NTVMHVDTVY DRPSNTTTEI WTLPEVLGER YSADKDVVVL TSGHTGGVAE
DIAYILKQMR RAIVVGERTE GGALDLQKLR IGQSNFFLTV PVSRLGPLG GGGQTWEGSG
VLPCVGTPAE QALEKALAIL

[REFERENCES]

1. Taniguchi R.T., *et al.* (2012) Proc. Natl. Acad. Sci. U.S.A. 109:7847-7852.
2. Zhou R., *et al.* (2012) J. Immunol. 188:1742-1750.
3. Jelcick A.S., *et al.* (2011) PLoS ONE 6:e21858-e21858.
4. Wisard J., *et al.* (2011) Invest. Ophthalmol. Vis. Sci. 52:5804-5811.